

AD-A107 628

LETTERMAN ARMY INST OF RESEARCH PRESIDIO OF SAN FRANC--ETC F/8 6/20
THE MUTAGENIC POTENTIAL OF: (E)-1,2,3,4-TETRAHYDRO-6-METHYL-1(2--ETC(U)
SEP 81 L J SAUERS, F R PULLIAM, J T FRUIN

UNCLASSIFIED

LAIR-109

NL

1 of 1
AD
A1076

END
DATE
FILMED
1 82
DTIC

AD A107628

LEVEL

2

INSTITUTE REPORT NO. 109

THE MUTAGENIC POTENTIAL OF:

(E)-1,2,3,4-tetrahydro-6-methyl-1-(2-methyl-1-oxo-2-butenyl)quinoline
1,2,3,4-tetrahydro-6-methyl-1-(3-methyl-1-oxo-2-butenyl)quinoline
50% DEET, 25% Dow Corning 200 fluid, in isopropanol .

LEONARD J. SAUERS/BA, SP5
FREDDICA R. PULLIAM/BS, SSG

JOHN T. FRUIN/DVM, PhD, LTC VC

TOXICOLOGY GROUP,
DIVISION OF RESEARCH SUPPORT

DTIC
ELECTE
NOV 20 1981


H

DISTRIBUTION STATEMENT A

Approved for public release;
Distribution Unlimited

SEPTEMBER 1981

Toxicology Series 20



LETTERMAN ARMY INSTITUTE OF RESEARCH PRESIDIO OF SAN FRANCISCO CALIFORNIA 94129

8 11 19 087

DTIC FILE COPY

Toxicology Series 20

Reproduction of this document in whole or in part is prohibited except with the permission of the Commander, Letterman Army Institute of Research, Presidio of San Francisco, California 94129. However, the Defense Technical Information Center is authorized to reproduce the document for United States Government purposes.

Destroy this report when it is no longer needed. Do not return it to the originator.

Citation of trade names in this report does not constitute an official endorsement or approval of the use of such items.

This material has been reviewed by Letterman Army Institute of Research and there is no objection to its presentation and/or publication. The opinions or assertions contained herein are the private views of the author(s) and are not to be construed as official or as reflecting the views of the Department of the Army or the Department of Defense. (AR 360-5)

John H. Marshall 23 Sept 87
(Signature and date)

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER LAIR Institute Report No. 109	2. GOVT ACCESSION NO. AD A107628	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) The Mutagenic Potential of: (E)-1,2,3,4-tetrahydro-6-methyl-1-(2-methyl-1-oxo-2-butenyl)quinoline; 1,2,3,4-tetrahydro-6-methyl-1-(3-methyl-1-oxo-2-butenyl)quinoline; 50% DEET, 25% Dow Corning 200 Fluid, in isopropanol.		5. TYPE OF REPORT & PERIOD COVERED Final June 1981 - September 1981
7. AUTHOR(s) Leonard J. Sauers, BA, SP5; Freddica R. Pulliam, BS, SSG; John T. Fruin DVM, PhD, LTC VC		6. PERFORMING ORG. REPORT NUMBER
9. PERFORMING ORGANIZATION NAME AND ADDRESS Toxicology Group, Div. of Research Support Letterman Army Institute of Research Presidio of San Francisco, CA 94129		6. CONTRACT OR GRANT NUMBER(s)
11. CONTROLLING OFFICE NAME AND ADDRESS U.S. Army Medical Research and Development Command Fort Detrick Frederick, MD 21701		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS Project 3M16770A871 Prevention of Military Disease Hazards WII 201
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		12. REPORT DATE September 1981
		13. NUMBER OF PAGES 36
		15. SECURITY CLASS. (of this report) UNCLASSIFIED
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) THIS DOCUMENT HAS BEEN APPROVED FOR PUBLIC RELEASE AND SALE: ITS DISTRIBUTION IS UNLIMITED.		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Mutagenicity, Toxicology, Ames Assay, (E)-1,2,3,4-tetrahydro-6-methyl-1-(2-methyl-1-oxo-2-butenyl)quinoline; 1,2,3,4-tetrahydro-6-methyl-1-(3-methyl-1-oxo-2-butenyl)quinoline; 50% DEET, 25% Dow Corning 200 Fluid, in isopropanol;		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) The mutagenic potential of (E)-1,2,3,4-tetrahydro-6-methyl-1-(2-methyl-1-oxo-2-butenyl)quinoline (CHR5); 1,2,3,4-tetrahydro-6-methyl-1-(3-methyl-1-oxo-2-butenyl)quinoline (CHR6); 50% DEET, 25% Dow Corning 200 Fluid, in isopropanol (CHF1); was assessed using the Ames/Salmonella Mammalian Microsome Mutagenicity Assay. Tester strains TA 98, TA 100, TA 1535, TA 1537, and TA 1538, were exposed to 1 ul/plate through 3.2×10^{-4} ul/plate doses of CHR5 and CHR6 and 0.1 ul/plate through 3.2×10^{-5} ul/plate doses of CHF1. No evidence of mutagenic activity was observed.		

DTIC
ELECTE
NOV 20 1981

H

DD FORM 1473

1 JAN 73

EDITION OF 1 NOV 65 IS OBSOLETE

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

ABSTRACT

The mutagenic potential of (E)-1,2,3,4-tetrahydro-6-methyl-1-(2-methyl-1-oxo-2-butenyl)quinoline(CHR5*); 1,2,3,4-tetrahydro-6-methyl-1-(3-methyl-1-oxo-2-butenyl)quinoline (CHR6*); 50% DEET, 25% Dow Corning 200 Fluid, in isopropanol (CHF1*) was assessed by using the Ames Salmonella/Mammalian Microsome Mutagenicity Assay. Tester strains TA 98, TA 100, TA 1535, TA 1537, and TA 1538 were exposed to 1 ul/plate through 3.2×10^{-4} ul/plate doses of CHR5 and CHR6 and 0.1 ul/plate through 3.2×10^{-5} ul/plate doses of CHF1. No evidence of mutagenic activity was observed.

* Code number for compound.

Accession For	
NTIS	CH43I
DTIC	IAB
Un-announced	
Justification	
By	
Distribution/	
Availability Codes	
Dist	Special
A	

PREFACE

AMES ASSAY REPORT:

<u>SUBSTANCE</u>	<u>CODE NO.</u>
(E)-1,2,3,4-tetrahydro-6-methyl-1-(2-methyl-1-oxo-1-butenyl)quinoline	CHR5
1,2,3,4-tetrahydro-6-methyl-1-(3-methyl-1-oxo-2-butenyl)quinoline	CHR6
50% DEET, 25% Dow Corning 200 Fluid in isopropanol	CHF1

TESTING FACILITY: Letterman Army Institute of Research
Presidio of San Francisco, CA 94129

SPONSOR: Division of Cutaneous Hazards
Letterman Army Institute of Research

PROJECT: More Effective Topical Repellents Against Disease Bearing Mosquitoes 3M62272A810

GLP STUDY NUMBER: 81017

STUDY DIRECTOR: LTC John T. Fruin D.V.M., PhD.

CO-PRINCIPAL INVESTIGATORS: SSG Freddica R. Pulliam, B.S.
SP5 Leonard J. Sauers, B.A.

RAW DATA: A copy of the final report, study protocol and retired SOPs will be maintained in the LAIR archives. Test substances were provided by sponsor. Chemical, analytical, stability, purity, etc. data are available from the sponsor.


PURPOSE: To determine the mutagenic potential of (E)-1,2,3,4-tetrahydro-6-methyl-1-(2-methyl-1-oxo-2-butenyl)quinoline; 1,2,3,4-tetrahydro-6-methyl-1-(3-methyl-1-oxo-2-butenyl)quinoline; 50% DEET, 25% Dow Corning 200 Fluid, in isopropanol, by using the Ames Salmonella/Mammalian Microsome Mutagenicity Test. Tester strains TA 98, TA 100, TA 1535, TA 1537, and TA 1538 were used.


ACKNOWLEDGMENTS


The authors wish to thank SP4 Thomas Kellner, BA; SP4 Larry Mullen, BS; and John Dacey for assistance in performing the research.

Signatures of Principal Scientists
Involved in the Study

We, the undersigned, believe the study, GLP number 81017, described in this report to be scientifically sound and the results and interpretation to be valid. The study was conducted to comply to the best of our ability with the Good Laboratory Practice Regulations outlined by the Environmental Protection Agency.


FREDDICA R. PULLIAM, BS Date
SSG
Co-Investigator


JOHN T. FRUIN, DVM, PhD Date
LTC, VC
Study Director


LEONARD J. SAUERS, BA Date
SP5
Co-Investigator



LETTERMAN 1950 CONCEPTS OF RESEARCH
PRESIDIO OF SAN FRANCISCO, CALIFORNIA 94129

REPLY TO
ATTENTION OF:

SGRD-ULZ-QA

22 July 1981

MEMORANDUM FOR RECORD

SUBJECT: Report of GLP Compliance

I hereby certify that in relation to IATR GLP study 81017 the following inspections were made:

1000 hr, 5 June 1981

1300 hr, 5 June 1981

Routine inspections with no adverse findings are reported quarterly, thus these inspections are also included in the July 1981 report to management and the Study Director.

JOHN C. JOHNSON
CPT, MS
Quality Assurance Officer

TABLE OF CONTENTS

Abstract.....	1
Preface.....	111
Acknowledgments.....	iv
Signatures of Principal Scientists.....	v
Report of Quality Assurance Unit.....	vi
Table of Contents.....	vii
BODY OF REPORT	
INTRODUCTION	
Rationale for using the Ames Assay.....	1
Description of Test, Rationale for strain selection.....	1
Description of Strains, History, Methods, and Data.....	2
METHODS	
Rationale for Dosage Levels and Response Tabulations....	3
Test Format.....	3
Statistical Analysis.....	4
RESULTS AND DISCUSSION.....	4
CONCLUSION.....	5
RECOMMENDATION.....	6
REFERENCES.....	7
APPENDIX (Tables 1 through 6).....	9
DISTRIBUTION LIST.....	25

Rationale for using the Ames Assay

The Ames Salmonella/Mammalian Microsome Mutagenicity Test is one of a standard bank of tests used by our laboratory for the assessment of the mutagenic potential of a test substance. It is a short-term screening assay for the prediction of potential mutagenic agents in mammals. It is inexpensive when compared to in vivo tests, yet is highly predictive and reliable in its ability to detect mutagenic activity and therefore carcinogenic probability (1). It relies on basic genetic principles and allows for the incorporation of a mammalian microsome enzyme system to increase sensitivity through enzymatically altering the test substance into an active metabolite. It has proven highly effective in assessing human risk (1).

Description of Test (Rationale for the selection of strains)

The test was developed by Bruce Ames, Ph.D. from the University of California-Berkeley. The test involves the use of several different genetically altered strains of Salmonella typhimurium, each with a specific mutation in the histidine operon (2). The test substance demonstrates mutagenic potential if it is able to revert the mutation in the bacterial histidine operon back to the wild type and thus reestablish prototrophic growth within the test strain. This reversion also can occur spontaneously due to a random mutational event. If, after adding a test substance, the number of revertants is significantly greater than the spontaneous reversion rate, then the test substance physically altered the locus involved in the operon's mutation and is able to induce point mutations and genetic damage (2).

In order to increase the sensitivity of the test system, two other mutations in the Salmonella are used (2). To insure a higher probability of uptake of test substance, the genome for the lipopolysaccharide layer (LP) is mutated and allows larger molecules to enter the bacteria. Each strain has another induced mutation which causes loss of excision repair mechanisms. Since many chemicals are not by themselves mutagenic but have to be activated by an enzymatic process, a mammalian microsome system is incorporated. These microsomal enzymes are obtained from livers of rats induced with Aroclor 1254; the enzymes allow for the expression of the metabolites in the mammalian system. This activated rat liver microsomal enzyme homogenate is termed S-9.

Description of Strains (History of the strains used, methods to monitor the integrity of the organisms, and data pertaining to current and historical controls and spontaneous reversion rates)

The test consists of using 99 strains of Salmonella typhimurium that are unable to grow without histidine because of a specific mutation in the histidine operon. This histidine requirement is verified by attempting to grow the tester strains on minimal glucose agar (MGA) plates, both with and without histidine. The dependence on this amino acid is shown when growth occurs only in its presence. The plasmids in strains TA 98 and TA 100 contain an ampicillin resistant R factor. Strains deficient in this plasmid demonstrate a zone of growth inhibition around an ampicillin impregnated disc. The alteration of the LP layer allows uptake by the Salmonella of larger molecules. If a crystal violet impregnated disc is placed onto a plate containing any one of the bacterial strains, a zone of growth inhibition will occur because the LP layer is altered. The absence of excision repair mechanisms can be determined by using ultraviolet (UV) light. These mechanisms function primarily by repairing photodimers between pyrimidine bases; exposure of bacteria to UV light will activate the formation of these dimers and cause cell lethality, since excision of these photodimers can not be made. The genetic mutation resulting in UV sensitivity also induces a dependence by the Salmonella to biotin. Therefore, this vitamin must be added. In order to prove that the bacteria are responsive to the mutation process, positive controls are run with known mutagens. If after exposure to the positive control substance, a larger number of revertants are obtained, then the bacteria are adequately responsive. Sterility controls are performed to determine the presence of contamination. Sterility of the test compound is also confirmed in each first dilution. Verification of the tester strains occurs spontaneously with the running of each assay. The value of the spontaneous reversion rate is obtained using the same inoculum of bacteria that is used in the assay (3).

Strains were obtained directly from Dr. Ames, University of California, Berkeley, and were maintained at -80 C in our laboratory. Before any substance was tested, quality controls were run on the bacterial strains to establish the validity of their special features and also to determine the spontaneous reversion rate (2). Records are maintained of all the data, to determine if deviations from the set trends have occurred.

We compared the spontaneous reversion values with our own historical values and those cited by Ames et al (2). Our conclusions are based on the spontaneous reversion rate compared to the experimentally induced rate of mutation. When operating effectively, these strains detect substances that cause base pair

mutations (TA 1535, TA 100) and frameshift mutations (TA 1537, TA 1538 and TA 98) (2).

METHODS (3)

Rationale for Dosage Levels and Dose Response Tabulations

To insure readable and reliable results, a sublethal concentration of the test substance had to be determined. This toxicity level was found by using MGA plates, various concentrations of the substance, and approximately 10^8 cells of TA 100 per plate, unless otherwise specified. Top agar containing trace amounts of histidine and biotin were placed on MGA plates. TA 100 is used because it is the most sensitive strain. Strain verification was confirmed on the bacteria, along with a determination of the spontaneous reversion rate. After incubation, the growth was observed on the plates. (The auxotrophic Salmonella will replicate a few times and potentially express a mutation. When the histidine and biotin supplies are exhausted, only those bacteria that reverted to the prototrophic phenotype will continue to reproduce and form macrocolonies; the remainder of the bacteria comprises the background lawn. The minimum toxic level is defined as the lowest serial dilution at which decreased macrocolony formation, below that of the spontaneous revertant rate, and an observable reduction in the density of the background lawn occurs.) A maximum dose of 1 mg/plate is used when no toxicity is observed. The densities were recorded as normal slight, and no growth.

Test Format

After we validated our bacterial strains and determined the optimal dosage of the test substance, we began the Ames Assay. In the actual experiment, 0.1ml of the particular strain of Salmonella (10^8 cells) and the specific dilutions of the test substance were added to 2 ml of molten top agar, which contained trace amounts of histidine and biotin. Since survival is better from cultures which have just passed the log phase, the Salmonella strains were used 16 hours (maximum) after initial inoculation into nutrient broth. The dose of the test substance spanned more than a 1000- fold, decreasing from the minimum toxic level by a dilution factor of 5. All the substances were tested with and without S-9 microsome fraction. The S-9 mixture which was previously titrated at an optimal strength was added to the molten top agar. After all the ingredients were added, the top agar was vortexed, then overlaid on minimum glucose agar plates. These plates contained 2% glucose and Vogel Bonner "E" Concentrate (4). The water used in this medium and all reagents came from a polymetric system. Plates were incubated, upside down in the dark at 37 C for 48 hours. Plates were prepared in triplicate and the average revertant counts were recorded. The corresponding number of revertants obtained was compared to the number of spontaneous

revertants; the conclusions were recorded statistically. A correlated dose response is considered necessary to declare a substance as a mutagen. Commoner (5), in his report, "Reliability of Bacterial Mutagenesis Techniques to Distinguish Carcinogenic and Non-Carcinogenic Chemical," and McCann et al (1) in their paper, "Detection of Carcinogens as Mutagen: Assay of over 300 Chemicals," have concurred on the test's ability to detect mutagenic potential.

Statistical Analysis

Quantitative evaluation was ascertained by two independent methods. Ames et al (2) assumed that a compound which caused twice the spontaneous reversion rate is mutagenic. Commoner (5), developed the MUTAR Ratio, which is stated in the following equation:

$$\text{MUTAR} = (E - C)/C_{AV}$$

Here, C is the number of spontaneous revertant colonies on control plates obtained on the same day and with the same treatment and strains. E is the number of revertants in response to the compound; C_{AV} is the number of spontaneous revertants on control plates calculated from historical records. The explanation of the results of this equation can be determined by the method of Commoner (5). This variation determines the probability of correctly classifying substances as carcinogens on the basis of their mutagenic activity. The E values were recorded by strain, with and without S-9. Values for C and C_{AV} were recorded separately.

We used the formula and logged all values for our permanent records.

RESULTS AND DISCUSSION

Throughout this report, all the test substances will be referred to by their respective code numbers:

<u>Substance</u>	<u>Code No</u>
(E)-1,2,3,4-tetrahydro-6-methyl-1-(2-methyl-1-oxo-2-butenyl)quinoline	CHR5
1,2,3,4-tetrahydro-6-methyl-1-(3-methyl-1-oxo-2-butenyl)quinoline	CHR6
50% DEET, 25% Dow Corning 200 Fluid, in Isopropanol	CHP1

On 3 June 1981, the Toxicity Level Determination was performed on the 3 test chemicals. For this experiment, all sterility, strain

verification, positive and negative controls were normal (Table 1). The plates containing the initial dilution showed slight growth for CHR5, normal growth for CHR6, and no growth for CHF1 (Tables 2A-2C). It was decided to use 0.1 ul/plate as the starting point for CHF1 and 1 ul/plate for the other test substances.

On 2 July 1981, The Ames Assay was run on the 3 test compounds. All sterility and strain verification controls were normal (Table 3). Unexpected results were observed in response to positive control chemical DMBA for all the strains. Expected results were seen in response to MNNG, AF, and BP, which validates our data since these other controls function through similar mechanisms. The negative controls were normal (Table 4).

No mutagenic activity was observed in response to test chemical CHR5 (Table 5A). One isolated incidence of mutagenic activity was seen for CHR6. This occurred at the 0.008 ul/plate dose for nonactivated TA 1537. No dose response was observed (Table 5B). A doubling of the spontaneous revertant rate was noticed in response to CHF1 at the 1.6×10^{-4} ul/plate level for nonactivated TA 1537 and nonactivated TA 1538, at the 0.02 ul/plate dose. No dose response was seen in either case (Table 5C).

The MUTAR values are listed in Tables 6A-6C. All calculations resulted in expected responses except for nonactivated TA 1538 at the 0.02 ul/plate dose level of CHF1 (Table 6C).

CONCLUSION

The results showed several isolated incidences of a doubling of the spontaneous reversion rate. It is in the opinion of the Ames Assay Laboratory at the University of California, Berkeley, that even though a doubling had occurred, one cannot declare mutagenicity unless an obvious dose response is seen (Maron D., Ames Assay Laboratory, University of California, Berkeley, 30 March 1981). Therefore on the basis of the Ames Test, compounds CHR5, CHR6, and CHF1 are not mutagenic at the levels tested.

RECOMMENDATION

We recommend that candidate insect repellents CHR5, CHR6, and CHF1 be tested further with other toxicological assays if efficacy tests show these compounds to be promising repellents.

REFERENCES

1. McCANN, J., E. CHOI, E. YAMASAKI, and B. N. AMES. Detection of carcinogens as mutagens in the Salmonella/microsome test: Assay of 300 chemicals. Proc Nat Acad Sci, USA 72:5135-5139, 1975
2. AMES, B. N., J. McCANN and E. YAMASAKI. Methods for detection carcinogens and mutagens with Salmonella/mammalian microsome mutagenicity test. Mutation Res 31: 347-364, 1975
3. LAIR SOP OP-STX-1, Ames Salmonella/mammalian microsome mutagenicity test, 1 March 1981
4. VOGEL, H. J. and D. M. BONNER. Acetylornithinase of E. coli: Partial purification and some properties. J Biol Chem 218: 97-106, 1956
5. COMMONER, B. Reliability of the bacterial mutagenesis techniques to distinguish carcinogenic and non-carcinogenic chemicals. EPA 600/1 76-022, 1976

LIST OF TABLES

	Date	Page
Table 1 Strain Verification for Toxicity Level Determination	3 June 1981	10
Table 2A Toxicity Level Determination	3 June 1981	11
Table 2B Toxicity Level Determination	3 June 1981	12
Table 2C Toxicity Level Determination	3 June 1981	13
Table 3 Strain Verification Control	2 July 1981	14
Table 4 Quality and Positive Controls	2 July 1981	15
Table 5A Salmonella/Microsome Assay Worksheet	2 July 1981	16
Table 5B Salmonella/Microsome Assay Worksheet	2 July 1981	18
Table 5C Salmonella/Microsome Assay Worksheet	2 July 1981	20
Table 6A Mutagenic Activity Ratio Worksheet	3 Aug 1981	22
Table 6B Mutagenic Activity Ratio Worksheet	3 Aug 1981	23
Table 6C Mutagenic Activity Ratio Worksheet	3 Aug 1981	24

APPENDIX

TABLE 1
STRAIN VERIFICATION FOR TOXICITY LEVEL DETERMINATION

Strains	Histidine Requirement	Ampicillin Resistance	UV	Sensitivity to Crystal Violet	Sterility Control	Response (1)
100	NG	G	NG	15.46 mm	NG	+
1537	NG	NG	NG	14.71 mm	NG	+
WT	G	NA	G	NA	NG	+

STERILITY CONTROL

His-Bio Mix Initial: NG End: NG MCA Plate: NG

Top Agar Initial: NG End: NG

Diluent: NG Nutrient Broth: NG

Test Compound (a) CHR5-113 (b) CHR6-NG (c) CHF1-NG (d) NA (e) NA

G = Growth NC = No Growth NT = Not Tested NA = Not Applicable WT = Wild Type

Spontaneous Revertants: TA 100, No S-9 - Average - 140 Positive Control - MNNG - 1612

(1) + = expected response - = unexpected response

Study Number: 81017 Date: 3 Jun 81 By: Sauers, Pulliam, Dacey, Mullen

TABLE 2A
TOXICITY LEVEL DETERMINATION

Substance assayed: CHR5 Substance dissolved in: EIOH
 Study Number: 81017 Date: 3 June 1981 Performed by: Sauers, Pulliam, Dacey, Mullen

TA 100 REVERTANT PLATE COUNT

Test Compound Concentration	Plate #1	Plate #2	Plate #3	Average	Background Lawn (1)
1 ul/plate	122	168	150	147	ST
0.1 ul/plate	145	167	188	167	NL
10 ⁻² ul/plate	132	121	167	140	NL
10 ⁻³ ul/plate	205	175	145	175	NL
10 ⁻⁴ ul/plate	137	147	129	138	NL
10 ⁻⁵ ul/plate	160	125	120	135	NL
10 ⁻⁶ ul/plate	150	180	168	166	NL
10 ⁻⁷ ul/plate	152	179	199	177	NL

(1) NC = No Growth ST = Slight Growth NL = Normal Lawn

TABLE 28

TOXICITY LEVEL DETERMINATION

Substance assayed: CHRG Substance dissolved in: ETOH
 Study Number: 81017 Date: 3 June 1981 Performed by: Sauers, Pulliam, Dacey, Mullen

TA 100 REVERTANT PLATE COUNT

Test Compound Concentration	Plate #1	Plate #2	Plate #3	Average	Background Lawn (1)
1 ul/plate	107	116	122	115	NL
0.1 ul/plate	88	97	111	99	NL
10 ⁻² ul/plate	111	90	104	102	NL
10 ⁻³ ul/plate	120	121	105	115	NL
10 ⁻⁴ ul/plate	131	144	115	130	NL
10 ⁻⁵ ul/plate	143	134	137	138	NL
10 ⁻⁶ ul/plate	134	162	135	144	NL
10 ⁻⁷ ul/plate	165	125	143	144	NL

(1) NG = No Growth ST = Slight Growth NL = Normal Lawn

TABLE 2C

TOXICITY LEVEL DETERMINATION

Substance assayed: Code CHF1 Substance dissolved in: ETOH
 Study Number: 81017 Date: 3 June 1981 Performed by: Sauers, Pulliam, Dacey, Mullen

TA 100 REVERTANT PLATE COUNT

Test Compound Concentration	Plate #1	Plate #2	Plate #3	Average	Background Lawn (1)
1 ul/plate	Toxic	Toxic	Toxic	Toxic	NG
0.1 ul/plate	85	86	110	94	NL
10 ⁻² ul/plate	136	119	102	119	NL
10 ⁻³ ul/plate	117	134	150	134	NL
10 ⁻⁴ ul/plate	123	119	158	133	NL
10 ⁻⁵ ul/plate	156	128	155	146	NL
10 ⁻⁶ ul/plate	144	158	115	139	NL
10 ⁻⁷ ul/plate	141	177	133	150	NL

(1) NG = No Growth ST = Slight Growth NL = Normal Lawn

TABLE 3

STRAIN VERIFICATION CONTROL

Strains	Histidine Requirement	Ampicillin Resistance	UV	Sensitivity to Crystal Violet	Sterility Control	
					Control	Response (1)
98	NG	G	NG	14.55 mm	NG	+
100	NG	G	NG	16.0 mm	NG	+
1535	NG	NT	NG	15.80 mm	NG	+
1537	NG	25.25 mm	NG	16.47 mm	NG	+
1538	NG	NT	NG	15.42 mm	NG	+
WT	G	NT	G	NT	NA	+

STERILITY CONTROL

His-Bio Mix Initial: NG End: NG Diluent: NG
 Top Agar Initial: NG End: NG MGA Plate: NT
 S-9 Mix Initial: NG End: NG Nutrient Broth: NT
 Test Compound (a) CHF1=NG (b) CHR5=NG (c) CHR6=NG (d) NA (e) NA (f) NA
 G = Growth NG = No Growth NT = Not Tested NA = Not Applicable WT = Wild Type
 Study Number: 81017 By: Dacey, Pulliam,
 Date: 2 July 1981 Kellner
 (1) + = expected response
 - = unexpected response

TABLE 4

POSITIVE CONTROL REVERTANT RATE AND SPONTANEOUS REVERTANT RATE

Compd.	Amount of Compd. Added	S-9 Added	98	100	Strain Number	
					1535	1537
AF	2 ug/plate	yes	(360,413,270) (348)	(243,426,215) (295)	(55,Tox,Tox) (55)	(55,Tox,Tox) (55)
BP	2 ug/plate	yes	(100,68,52) (73)	(446,267,271) (328)	(53,26,22) (34)	(Tox,Tox,Tox) (Tox)
DMBA	20 ug/plate	yes	(19,38,20) (26)	(170,131,166) (156)	(10,18,15) (14)	(15,9,13) (12)
MMNG	2 ug/plate	no		(985,1041,1050) (1025)		
	20 ug/plate	no			(580,Tox,Tox) (580)	
<u>Spontaneous Revertant Test</u>						
	before		(30,14,24)	(122,101,122)	(13,16,15)	(12,11,7)
	after	yes	(12,20,8) (18)	(77,69,87) (96)	(14,8,5) (12)	(3,11,4) (8)
	before		(12,23,9)	(89,70,85)	(15,7,13)	(9,2,5)
	after	no	(9,18,7) (11)	(117,71,75) (85)	(11,17,15) (13)	(3,3,0) (4)
						(7,9,16)
						(8,NG*,NG*) (15)
						(11,1,NG*) (9)

Study Number: 81017

Date: 2 July 1981 By: Dacey, Kellner, Pulliam

* No background lawn

TABLE 5A
NUMBER OF REVERTANTS/PLATE

Compd.	Amount of Compd. Added	S-9 Added	Strain Number			
			98	100	1535	1537
CHR5	1 ul/plate	no	(10,16,7) (11)	(54,60,58) (57)	(12,9,12) (11)	(7,4,6) (6)
		yes	(17,28,16) (20)	(109,119,92) (107)	(19,17,9) (15)	(3,8,0) (4)
CHR5	0.2 ul/plate	no	(2,8,4) (5)	(63,83,82) (76)	(7,13,7) (9)	(2,5,3) (3)
		yes	(20,10,16) (15)	(59,68,79) (69)	(12,13,20) (15)	(7,6,3) (5)
CHR5	0.04 ul/plate	no	(2,14,7) (8)	(87,64,68) (73)	(22,8,9) (13)	(3,4,9) (5)
		yes	(27,24,24) (25)	(79,72,71) (74)	(16,13,21) (17)	(9,3,3) (5)

-continued

TABLE 5A, concluded
NUMBER OF REVERTANTS/PLATE

Compd.	Amount of Compd. Added	S-9 Added	NUMBER OF REVERTANTS/PLATE			
			98	100	Strain Number 1535 1537	1538
CHR5	0.008 ul/plate	no	(20,8,12) (13)	(98,60,71) (76)	(15,8,9) (11)	(6,3,3) (4)
		yes	(22,12,3) (12)	(84,106,95) (95)	(13,23,12) (16)	(8,8,6) (7)
CHR5	0.0016 ul/plate	no	(11,15,11) (12)	(62,87,78) (76)	(19,12,13) (15)	(9,3,6) (6)
		yes	(24,14,17) (18)	(66,57,55) (59)	(12,13,4) (10)	(5,8,4) (6)
CHR5	0.00032 ul/plate	no	(9,12,16) (12)	(81,79,61) (74)	(10,12,14) (12)	(3,6,2) (4)
		yes	(36,27,20) (28)	(110,80,94) (95)	(11,13,11) (12)	(7,3,6) (5)

Study Number: 81017

Date: 2 July 1981

By: Dacey, Kellner, Pulliam

TABLE 5B
NUMBER OF REVERTANTS/PLATE

Compd.	Amount of Compd. Added	S-9 Added	Strain Number			
			98	100	1535	1538
CHR6	1 ul/plate	no	(Tox,Tox,Tox) (Tox)	(Tox,Tox,Tox) (Tox)	(Tox,Tox,Tox) (Tox)	(Tox,Tox,Tox) (Tox)
		yes	(Tox,Tox,23) (Tox)	(55,Tox,66) (61)	(Tox,Tox,Tox) (Tox)	(Tox,Tox,Tox) (Tox)
CHR6	0.2 ul/plate	no	(15,17,11) (14)	(88,26,72) (62)	(6,12,6) (8)	(9,8,9) (9)
		yes	(23,10,21) (18)	(105,126,114) (115)	(26,19,13) (21)	(1,3,1) (13)
CHR6	0.04 ul/plate	no	(14,10,10) (15)	(78,66,88) (77)	(6,15,6) (9)	(5,4,2) (5)
		yes	(22,30,20) (24)	(107,106,85) (99)	(14,15,13) (14)	(1,19,2a) (12)
CHR6	0.008 ul/plate	no	(21,13,12) (15)	(76,27,92) (85)	(6,22,5) (11)	(10,3,11) (10)
		yes	(20,18,12) (17)	(67,72,82) (74)	(12,8,9) (10)	(19,16,7) (14)

-continued

Study Number: 81017

Date: 2 July 1981

By: Dacey, Kellner, Pulliam

TABLE 5B, concluded
NUMBER OF REVERTANTS/PLATE

Compd.	Amount of Compd. Added	S-9 Added	Strain Number			
			98	100	1535	1537
CHR6	0.0016 ul/plate	no	(5,14,9) (9)	(57,73,70) (67)	(6,15,12) (11)	(1,17,6) (8)
		yes	(19,13,13) (15)	(78,78,61) (72)	(13,9,8) (10)	(5,5,6) (5)
CHR6	0.00032 ul/pl.	no	(5,19,10) (11)	(81,56,68) (68)	(4,12,10) (9)	(2,6,3) (4)
		yes	(11,10,20) (14)	(52,81,65) (66)	(8,6,6) (7)	(4,8,2) (5)
						(10,12,10) (11)
						(10,19,14) (14)

Study Number: 81017 Date: 2 July 1981 By: Dacey, Kellner, Pulliam

TABLE 5C
NUMBER OF REVERTANTS/PLATE

Compd.	Amount of Compd. Added	S-9 Added	Strain Number			
			98	100	1535	1537
CHF1	0.1 ul/plate	no	(21,16,12) (16)	(141,106,74) (107)	(15,10,17) (14)	(8,5,7) (7)
					(38,13,4) (18)	
		yes	(32,27,25) (28)	(150,123,110) (128)	(22,22,18) (21)	(8,3,16) (9)
					(11,13,33) (19)	
CHF1	0.02 ul/plate	no	(21,23,20) (21)	(110,111,113) (111)	(12,30,18) (20)	(3,4,5) (4)
		yes	(10,20,8) (13)	(70,70,88) (76)	(16,13,12) (14)	(10,6,13) (10)
					(5,11,6) (7)	
CHF1	0.004 ul/plate	no	(3,4,5) (4)	(73,64,66) (63)	(9,2,7) (3)	(3,4,4) (3)
		yes	(17,21,17) (18)	(95,81,64) (80)	(16,12,11) (13)	(5,2,5) (4)
					(20,23,16) (20)	
CHF1	0.0008 ul/pl.	no	(5,12,4) (7)	(64,56,65) (62)	(9,8,11) (9)	(9,10,6) (3)
		yes	(27,36,24) (29)	(99,93,86) (93)	(16,18,21) (18)	(9,4,4) (6)
					(14,14,17) (15)	

-continued

Study Number: 81017

Date: 2 July 1981

By: Dacey, Kellner, Pulliam

TABLE 5C, concluded
NUMBER OF REVERTANTS/PLATE

Compd.	Amount of Compd. Added	S-9 Added	98	100	Strain Number		
					1535	1537	1538
CHF1	0.00016 ul/pl.	no	(6,15,19) (13)	(98,93,93) (95)	(14,22,11) (16)	(11,12,9) (11)	(8,15,17) (13)
		yes	(16,18,25) (20)	(103,101,126) (110)	(22,12,18) (17)	(9,8,2) (6)	(19,20,18) (19)
CHF1	0.000032 ul/pl.	no	(24,22,14) (20)	(123,111,102) (112)	(18,15,8) (14)	(9,3,6) (6)	(13,7,11) (10)
		yes	(23,26,20) (23)	(81,66,120) (89)	(12,21,19) (17)	(6,3,6) (5)	(12,4,0) (5)

Study Number: 81017 Date: 2 July 1981 By: Dacey, Kellner, Pulliam

TABLE 6A
MUTAGENIC ACTIVITY RATIO

Substance Assayed: CHR5 Dissolved in: ETOH

Study Number: 81017 Date: 8 Aug 81 By: Sauers

Concentration	Strain	MUTAR (act)	MUTAR	Concentration	Strain	MUTAR (act)	MUTAR
1 ul/plate	TA 98	0.09	*	0.008 ul/pl.	TA 1535	0.36	*
0.2 ul/plate	TA 98	*	*	0.0016 ul/pl.	TA 1535	*	0.13
0.04 ul/plate	TA 98	0.30	*	0.00032 ul/pl.	TA 1535	*	*
0.008 ul/plate	TA 98	*	0.11				
0.0016 ul/pl.	TA 98	*	0.05	1 ul/plate	TA 1537	*	0.33
0.00032 ul/pl.	TA 98	0.43	0.05	0.2 ul/plate	TA 1537	*	*
				0.04 ul/plate	TA 1537	*	0.17
1 ul/plate	TA 100	0.1	*	0.008 ul/plate	TA 1537	*	*
0.2 ul/plate	TA 100	*	*	0.0016 ul/pl.	TA 1537	*	0.33
0.04 ul/plate	TA 100	*	*	0.00032 ul/pl.	TA 1537	*	*
0.008 ul/plate	TA 100	*	*				
0.0016 ul/pl.	TA 100	*	*	1 ul/plate	TA 1538	0.5	*
0.00032 ul/pl.	TA 100	*	*	0.2 ul/plate	TA 1538	0.17	*
				0.04 ul/plate	TA 1538	0.33	0.23
1 ul/plate	TA 1535	0.27	*	0.008 ul/pl.	TA 1538	*	*
0.2 ul/plate	TA 1535	0.27	*	0.0016 ul/pl.	TA 1538	*	*
0.04 ul/plate	TA 1535	0.45	*	0.00032 ul/pl.	TA 1538	*	*

(act): S-9 fraction was added

* : calculated value resulted in a negative MUTAR, or a zero MUTAR

TABLE 6B
MUTAGENIC ACTIVITY RATIO

Substance Assayed: CHR6 Dissolved in: ETOH

Study Number: 81017 Date: 3 August 81 By: Sauers

Concentration	Strain	MUTAR (act)	MUTAR	Concentration	Strain	MUTAR (act)	MUTAR
1 ul/plate	TA 98	*	*	0.008 ul/plate	TA 1535	*	*
0.2 ul/plate	TA 98	*	0.16	0.0016 ul/pl.	TA 1535	*	*
0.04 ul/plate	TA 98	0.26	0.21	0.00032 ul/pl.	TA 1535	*	*
0.008 ul/plate	TA 98	*	0.21				
0.0016 ul/pl.	TA 98	*	*	1 ul/plate	TA 1537	*	*
0.00032 ul/pl.	TA 98	*	*	0.2 ul/plate	TA 1537	*	*
				0.04 ul/plate	TA 1537	*	*
1 ul/plate	TA 100	*	*	0.008 ul/plate	TA 1537	*	0.83
0.2 ul/plate	TA 100	0.18	*	0.0016 ul/pl.	TA 1537	*	0.67
0.04 ul/plate	TA 100	0.03	*	0.00032 ul/pl.	TA 1537	*	*
0.008 ul/plate	TA 100	*	*				
0.0016 ul/pl.	TA 100	*	*	1 ul/plate	TA 1538	*	*
0.00032 ul/pl.	TA 100	*	*	0.2 ul/plate	TA 1538	*	*
				0.04 ul/plate	TA 1538	0.17	*
1 ul/plate	TA 1535	*	*	0.08 ul/plate	TA 1538	*	0.08
0.2 ul/plate	TA 1535	0.8	*	0.0016 ul/pl.	TA 1538	0.28	*
0.04 ul/plate	TA 1535	0.18	*	0.00032 ul/pl.	TA 1538	*	0.15

(act): S-9 fraction was added

* : calculated value resulted in a negative MUTAR, or a zero MUTAR

TABLE 6C
MUTAGENIC ACTIVITY RATIO

Substance Assayed: CHF1 Dissolved in: ETOH
Study Number: 81017 Date: 3 August 1981 By: Sauers

Concentration	Strain	MUTAR (act)	MUTAR	Concentration	Strain	MUTAR (act)	MUTAR
0.1 ul/plate	TA 98	0.43	0.26	0.0008 ul/pl.	TA 1535	0.54	*
0.02 ul/plate	TA 98	*	0.53	0.00016 ul/pl.	TA 1535	0.45	0.2
0.004 ul/plate	TA 98	*	*	0.000032 ul/pl.	TA 1535	0.45	0.07
0.0008 ul/pl.	TA 98	0.47	*				
0.00016 ul/pl.	TA 98	0.09	0.11	0.1 ul/plate	TA 1537	0.15	0.5
0.000032 ul/pl.	TA 98	0.21	0.47	0.02 ul/plate	TA 1537	0.29	*
				0.004 ul/plate	TA 1537	*	*
0.1 ul/plate	TA 100	0.3	0.24	0.0008 ul/pl.	TA 1537	*	1.17
0.02 ul/plate	TA 100	*	0.28	0.00016 ul/pl.	TA 1537	*	0.33
0.004 ul/plate	TA 100	*	*	0.000032 ul/pl.	TA 1537	*	0.33
0.0008 ul/pl.	TA 100	*	*				
0.00016 ul/pl.	TA 100	0.13	0.11	0.1 ul/plate	TA 1538	0.22	0.68
0.000032 ul/pl.	TA 100	*	0.29	0.02 ul/plate	TA 1538	*	1.59
				0.004 ul/plate	TA 1538	0.28	0.08
0.1 ul/plate	TA 1535	0.8	0.07	0.0008 ul/pl.	TA 1538	*	*
0.02 ul/plate	TA 1535	0.18	0.47	0.00016 ul/pl.	TA 1538	0.22	0.3
0.004 ul/plate	TA 1535	0.09	*	0.000032 ul/pl.	TA 1538	*	0.08

(act): S-9 fraction was added

* : calculated value resulted in a negative MUTAR, or a zero MUTAR

OFFICIAL DISTRIBUTION LIST

Comander
US Army Medical Research and Development Command
ATTN: SGRD-SI/ Mrs. Madigan
Fort Detrick, Frederick MD 21701

Defense Technical Information Center
ATTN: DTIC-DDA (12 copies)
Cameron Station
Alexandria VA 22314

Director of Defense Research and Engineering
ATTN: Assistant Director, Environmental and
Life Sciences
Washington DC 20301

The Surgeon General
ATTN: DASG-TLO
Washington DC 20314

HQ DA (DASG-ZXA)
WASH DC 20310

Superintendent
Academy of Health Sciences
ATTN: AHS-COM
Fort Sam Houston TX 78234

Assistant Dean
Institute and Research Support
Uniformed Services University of Health Sciences
6917 Arlington Road
Bethesda MD 20014

Commander
US Army Environmental Hygiene Agency
Aberdeen Proving Ground MD 21070

US Army Research Office
ATTN: Chemical and Biological Sciences Division
P.O. Box 1221
Research Triangle Park NC 27709
Biological Sciences Division
Office of Naval Research
Arlington VA 22217

Director of Life Sciences
USAF Office of Scientific Research (AFSC)
Bolling AFB
Washington DC 20332

Director
Walter Reed Army Institute of Research
Washington DC 20012

Commander
US Army Medical Research Institute of Infectious
Diseases
Fort Detrick, Frederick MD 21701

Commander
US Army Research Institute of Environmental
Medicine
Natick MA 01760

Commander
US Army Institute of Surgical Research
Brooke Army Medical Center
Fort Sam Houston TX 78234

Commander
US Army Institute of Dental Research
Washington DC 20012

Commander
US Army Medical Bioengineering
Research and Development Laboratory
Fort Detrick, Frederick MD 21701

Commander
US Army Aeromedical Research Laboratory
Fort Rucker AL 36362

Commander
US Army Biomedical Laboratory
Aberdeen Proving Ground
Edgewood Arsenal MD 21010

Commander
Naval Medical Research Institute
National Naval Medical Center
Bethesda MD 20014

Commander
USAF School of Aerospace Medicine
Aerospace Medical Division
Brooks Air Force Base TX 78235

**DAT
FILM**